


# Optimized PTFE Filmcast Liner Technology (O-PTFE)


# Flexibility Meets Durability.

Heat Shrink



	APPLICATION	BENEFIT	
Durability & Abrasion Resistance	Minimizes delivery forces while maintaining flexibility	Lower delivery forces and particle generation spec. relevant for stent retriever, flow diverter delivery	
High-Pressure Resistance	Accommodates the delivery of embolization agents	Optimized for aspiration and embolization catheters	
Ultra-Thin Walls	Maximizes flow rates and facilitates device interface	Minimized wall thickness for greater ID and distal flexibility	

#### **DURABILITY AND LUBRICITY**




Bench Testing in anatomical model measuring the retraction force of a stent retriever being pulled through the ID of a PTFE lined catheter.

## **RESULTS**

- Durability and Lubricity comparable to extruded PTFE
- Up to 30% improvement in durability and lubricity vs conventional filmcast.

#### **FLEXIBILITY**




Tensile testing performed to generate stress/strain curves for legacy and O-PTFE to generate tensile modulus to determine impact on material stiffness.

#### **RESULTS**

No trade-off in Flexibility vs. Conventional PTFE

#### **CONVENTIONAL PTFE VS O-PTFE BURST PRESSURE**



Liner only was tested to burst using water to pressurize the ID.

# **RESULTS**

- Up to 25% improvement over legacy PTFE
- Specifically suitable for aspiration and embolization catheters.

### **SPECIFICATIONS**

Mandrel	Inside	ID	Nominal	Wall Tolerance	Cut
	Diameter (ID)	Tolerance	Wall Thickness	(average)	Length
Annealed Stainless Steel/ Silver Plated Copper	0.008" - 0.100"	± 0.0002" - ± 0.0005"	0.0004" - 0.003"	± 0.0001" - ± 0.0005"	85" Max./Spooled

