

Oscor® Inc. – an Integer Company

# Destino™


## Steerable Guiding Sheath

The Destino™ is a next generation bi-directional steerable guiding sheath for the introduction and placement of diagnostic and therapeutic devices into the human vasculature, including but not limited to intracardiac, transseptal, and other placements.

Featuring a next generation hemostatic valve and braided FlexSteer™ Shaft Technology, the Destino™ is the only available steerable sheath in a size range of 8.5F, 10F, and 12F with true bi-directional tip deflections of 180°.








# **SHAFT DIAMETERS** 8.5 F 10 E 14.1 F 15.9 F 12.1 F 4.70 mm/ 4.10 mm / 5.30 mm / 0.185" 0.158" 0.208" TIP TIP 2.80 mm , 0.112" C: 17mm **CURVE CONFIGURATIONS** 8.5F 50 mm 61 cm 22 mm Integer®

## **ORDERING INFORMATION**





#### Destino™ Kit:

- (1) Steerable Guiding Sheath
- (1) Straight or 50  $^{\circ}$  Curved Dilator
- (1) Guidewire, OD: 0.035", L: 180 cm

#### Destino<sup>™</sup> Set:

- (1) Steerable Guiding Sheath
- (1) Straight or 50 ° Curved Dilator



		<b>A1</b>	A2	<b>A3</b>	B1	B2	<b>A4</b>	С
Model # Kit (One per box)	Model # Set (One per box)	Sheath Inner Diameter	Usable Sheath Length	Overall Sheath Length	Usable Dilator Length	Overall Dilator Length	Sheath Curve Size	Compatible Transseptal Needle
Straight Dilator								
DS857117GW	DS857117	8.5 F	71 cm	91 cm	94 cm	96 cm	17 mm	98 cm
DS857122GW	DS857122	8.5 F	71 cm	91 cm	94 cm	96 cm	22 mm	98 cm
DS857150GW	DS857150	8.5 F	71 cm	91 cm	94 cm	96 cm	50 mm	98 cm
DS106722GW	DS106722	• 10 F	67 cm	87 cm	90 cm	92 cm	22 mm	98 cm
DS126750GW	DS126750	<b>12</b> F	67 cm	87 cm	90 cm	92 cm	50 mm	98 cm
50° Curved Dilato	r							
DC857117GW	DC857117	8.5 F	71 cm	91 cm	94 cm	96 cm	17 mm	98 cm
DC857122GW	DC857122	8.5 F	71 cm	91 cm	94 cm	96 cm	22 mm	98 cm
DC857150GW	DC857150	8.5 F	71 cm	91 cm	94 cm	96 cm	50 mm	98 cm
DC106722GW	DC106722	<b>10</b> F	67 cm	87 cm	90 cm	92 cm	22 mm	98 cm
DC126750GW	DC126750	• 12 F	67 cm	87 cm	90 cm	92 cm	50 mm	98 cm

Guiding sheaths are supplied sterile. For single use only. Standard lengths are noted above. Other lengths available upon request.

**Caution:** Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.



P: +1 727-937-2511 sales@oscor.com | **www.oscor.com**