

Oscor® Inc. – an Integer Company

Adaptors & Extensions

Bifurcated, Unipolar, and Bipolar Implantable Lead Adaptors, Extensions and Plugs

Adaptors link existing leads to new active implantable medical devices. We select the best material and utilize leading industry processes to ensure optimal performance and reliability between new and existing technologies. In addition to our custom development programs, we offer non-proprietary adaptors to speed development timelines. Our standard configurations are sold through the major OEM companies as additions to product lines in the markets of cardiac rhythm, neuromodulation and hearing restoration.

BIFURCATED ADAPTORS

iLINK-BLV*



*Only available in USA

BIPOLAR ADAPTORS

BLV/BIS



iLINK-BIS*



*Only available in USA

BIS/BLV



2XBIS/BIS



B/IS



BIPOLAR EXTENSIONS

BIS/BIS



UNIPOLAR ADAPTORS

C/IS Adaptor



M/IS Adaptor



PLUGS

VV



LV-1



IS4/DF4





ORDERING INFORMATION

Model #	Receptacle	Connector	Coil	Insulation	Length
ADAPTORS	BIFURCATED	BIPOLAR			
iLINK-BLV-10*	One Bipolar LV-1 One Bipolar IS-1	IS-1	Quadrifilar	Silicone	10 cm
iLINK-BIS-10*	Two Bipolar IS-1	IS-1	Quadrifilar	Silicone	10 cm
2XBIS/BIS-17	Two Bipolar IS-1	IS-1	Quadrifilar	Silicone	17 cm
B/IS-15 SS2	Two Unipolar 5mm	IS-1	Quadrifilar	Silicone	15 cm
ADAPTORS	BIPOLAR	BIPOLAR			
BLV/BIS-10	LV-1	IS-1	Quadrifilar	Silicone	10 cm
BLV/BIS-4403	LV-1	IS-1	Quadrifilar	Silicone	17 cm
BLV/BIS-40	LV-1	IS-1	Quadrifilar	Silicone	40 cm
BIS/BLV-4402	IS-1	LV-1	Quadrifilar	Silicone	17 cm
EXTENSIONS	BIPOLAR	BIPOLAR			
BIS/BIS-17	IS-1	IS-1	Quadrifilar	Silicone	17 cm
BIS/BIS-40	IS-1	IS-1	Quadrifilar	Silicone	40 cm
ADAPTORS	UNIPOLAR	UNIPOLAR			
C/IS-10	6 mm	IS-1	Quadrifilar	Silicone	10 cm
C/IS-10 SS2	6 mm	IS-1	Quadrifilar	Silicone	10 cm
M/IS-10	5 mm	IS-1	Quadrifilar	Silicone	10 cm
PLUGS		CONNECTOR			
VV-PLUG	N/A	IS-1	N/A	Silicone	3.5 cm
LV-1 PLUG	N/A	LV-1	N/A	Silicone	3.4 cm
IS4/DF4 PLUG	N/A	IS4/DF4	N/A	Silicone	4 cm

Standard lengths are noted above. Other lengths available upon request. When ordering please add the usable length to model.
*Only available in USA. Products not available in all regions.

Caution: Federal Law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the Instructions for Use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.

